

— REUSABLE METAL PEN · MANUAL PRESS-INJECT

ES04

A reusable, manually driven injection pen in a premium anodized metal housing — one pen, up to 3 years or 3,000 cycles. Rotate the external dose ring, press to inject, and reuse. Fully customizable, one supply chain: pen, cartridge, needle and carton.




75 clicks Max dose per injection	0.01 mL Dose step (1 click)	3.0 mL Cartridge · ISO 11608-2	Metal Anodized housing	3,000 cycles Reusable · 3 years	1,000 pcs MOQ
--	---------------------------------------	--	----------------------------------	---	-------------------------

— Specifications

Injection method	Manual — external rotation to set dose, then press the actuator
Dose-setting method	External rotation (turn the external dose ring)
Available doses	60 clicks (0.60 mL) or 75 clicks (0.75 mL)
Dose step	0.01 mL per click
Built-in container	3.0 mL cartridge, ISO 11608-2 compliant
Usage type	Reusable — lifespan of 3 years or 3,000 cycles
Compatible needle	Sterile pen needle, ISO 11608-2 compliant
Body material & colour	Anodized aluminum; champagne gold, blue, green, maroon, magenta

DIMENSIONS



WHY CUSTOMERS CHOOSE IT

Metal housing Rigid, corrosion-resistant anodized aluminum that holds up over 3 years or 3,000 cycles.	Stable, accurate delivery The ergonomic plunger maintains consistent injection force with clear dosage feedback.
Wide cartridge compatibility Standard 3.0 mL ISO 11608-2 cartridges across insulin, GLP-1, interleukin and IgE biologics.	Reusable by design One pen serves for up to 3 years or 3,000 cycles, cutting cost-per-dose and waste.

THE DEVICE

Six anodized finishes	Single pen
Cap · body · cartridge	Profile detail

Waytop is a wholesale supplier of injection pens and related components. Regulatory market entry, local certifications and compliance are the buyer's responsibility. Waytop does not sell into markets where a given model is subject to an existing authorized distributor or patent. Intended applications, dose ranges and all clinical statements are determined by the final registered labeling for each market.